




1	bis Fert. Nr. 4812 22 0473/44 für GT 1900 22 0474/44 für GT 2500 22 0475/47 für GT 3200 22 0476/48 für GT 4000 22 0477/49 für GT 5000 ab Fert. Nr. 4901 22 0496/44 für GT 1900/GT 2000 22 0497/44 für GT 2500 22 0498/47 für GT 3200 22 0499/47 für GT 4000 22 0500/49 für GT 5000	10 5232/20 ab Fert.Nr. 5104 für GT 3200/4000 10 3650/21 bis Fert.-Nr. 5103 10 5230/20 ab Fert.Nr. 5104 für GT 5000 10 3843/21 bis Fert.-Nr. 5103 10 5234/20 ab Fert.Nr. 5104 01 5015/02 Gleitkörper im Scharnier 12 01 3485/05/25 13 10 3653/16 14 GT 1900 GT 2000 GT 2500 07 2131/18 GT 3200 GT 4000 07 2283/18 Aspera 07 2131/18 Bosch GT 5000 07 2327/18 15 Relais GT 1900 GT 2000 GT 2500 07 2130/18 bis Fert.-Nr. 5601 07 3899/09 ab Fert.-Nr. 5602 GT 3200 GT 4000 07 2093/18 Aspera 07 3043/16 Bosch 07 3044/18 Relaiskomb. Necchi GT 5000 07 2328/18 Aspera 07 2533/18 Relaiskomb. Necchi 16 04 4108/03	17 für GT 3200/4000 und GT 5000 04 4018/15 Metall 10 7995/03 Kunststoff 13 0455/28 für GT 3200/ GT 4000/5000 10 3530/05 Kappe 19 10 3602/02 für GT 1900/2000/2500 04 3338/07 für GT 3200/4000/5000 10 4548/08 für GT 5000 (Necchi) 20 04 1994/02 für GT 1900/2000/2500 und GT 5000 21 01 3486/03 22 13 0664/45 für GT 1900/2000 und GT 2500 04 2095/06 Platte 13 0307/45 GT 3200/GT 4000 13 0870/49 für GT 5000 13 0311/47 für GT 5000 (Necchi) 23 12 3696/18 24 Saugrohr 10 3646/26 für GT 1900/2000/2500 10 3651/29 für GT 3200/4000 10 3845/25 für GT 5000 10 4720/10 für GT 5000 (Necchi) 10 4767/11 für GT 3200/ 4000/5000 Saugrohr-Verlängerung (Necchi) 10 4768/06 Servicerohr für GT 5000 (Necchi)	25 07 3651/27 26 10 3847/12 bis Fert. Nr. 4901 10 4001/12 ab Fert. Nr. 4902 10 4796/14 ab Fert.Nr. 5209 27 01 3551/04/25 Blattfeder 28 04 4109/03/2 04 4020/04 mit Schlitz 29 Kontrolleuchte 04 3391/11 rot 04 3393/11 grün 04 3395/11 gelb 30 04 4110/12 Schalter 31 01 3552/05/2 bis Fert. Nr. 4812 01 3677/05/2 ab Fert. Nr. 4901 32 Schloß 04 4111/18 bis Fert. Nr. 4812 04 4341/15 ab Fert. Nr. 4901 04 5523/11 ab Fert.Nr. 5104 Stift für Befestigung Schloß 01 3624/04/25 bis Fert. Nr. 4812 01 3678/03/10 ab Fert. Nr. 4901 33 45 9007/99 kubisch 19 3497/21 konisch 34 19 1411/04 für GT 1900/2000 19 1413/05 für GT 2500 19 1415/05 für GT 3200 19 1417/06 für GT 4000 19 1773/06 für GT 5000	50 8671/00 Bed.-Anleitung 50 8202/02 Gefrierfibel
---	---	--	---	---	---